

MATERIAL SAFETY DATA SHEET

DELASCO INSTRUMENT LUBRICANT

Delasco Dermatologic Lab & Supply, Inc.
608 13th Ave., Council Bluffs, IA 51501 1-712-323-3269
www.delasco.com questions@delasco.com
Date prepared: December 2010

The statements herein are offered for informational purposes only and are intended to be followed only by persons having related technical skills and at their own discretion and risk. Since conditions and manner of use are outside our control, we make no warranties, express or implied, and assume no liability in connection with any use of this information.

Product/Chemical Name: DELASCO INSTRUMENT LUBRICANT

Chemical Description: Polymeric Lubricant

Physical Description: White, milk like liquid with slight oleic odor. Emulsifies in water. Specific gravity = 1.00.

Special Hazards: Low hazard. Not a D.O.T. hazardous material as defined by 49 CFR parts 171-173.

Unusual Fire/Explosion Hazards: None

Reactivity: Stable

Hazardous Decomposition Products: None known.

Spill or Leak Procedures

Action: Take up with absorbent material and place into disposal container.

Disposal Method: To be performed in compliance with all current local, state and federal regulations.

Overexposure Effects: Avoid prolonged or repeated skin contact.

Potential Carcinogen: Not available.

First Aid Procedures

Skin: First aid is typically not required. However, grossly contaminated clothing should be removed.

Eyes: Check for and remove any contact lenses. Wash eyes immediately with large amounts of water, occasionally lifting upper and lower eyelids, until no evidence of the material remains (at least 15-20 minutes). Get medical attention if irritation occurs.

Ingestion: First aid is typically not required. Not expected to cause serious harm.

Inhalation: If irritation occurs, move victim to fresh air and seek medical attention.

Storage Recommendations: Store at room temperature. Keep container tightly closed when not in use.

Special Protection Information

Protective eye goggles and rubber gloves recommended.

Special Handling Precautions

Handle with care. Avoid unnecessary contact. Wash thoroughly after handling.