

Electrosurgery Unit Comparison Chart						
	Delasco Electricator	Aaron 900	Aaron 950	Hyfrecator 2000	Ellman Radiowave	Ellman Surgitron EMC
Dimensions	8.5"(h) x 6.5"(w) x 3.0"(d)	6.3"(h) x 9"(w) x 4.5"(d)	7"(h) x 10.3"(w) x 4.5"(d)	7.5"(h) x 8.75"(w) x 4"(d)	5"(h) x 8"(w) x 5"(d)	6.25"(h) x 8"(w) x 9"(d)
Weight	6 lbs	5 lbs	7 lbs	6 lbs	7 lbs	9.25 lbs
Max Output	1-25 Watts (low) 1-45 Watts (high)	30 Watts	35 Watts	1-20 Watts (low) 1-35 Watts (high)	90 Watts	140 Watts
Cut	N/A	N/A	60 Watts	N/A	90 Watts	140 Watts
Coag	45 Watts	30 Watts	60 Watts	35 Watts	90 Watts	140 Watts
Blend	N/A	N/A	60 Watts	N/A	90 Watts	140 Watts
Bipolar Output	45 Watts	30 Watts	30 Watts	35 Watts	90 Watts	140 Watts
Frequency	350 - 490 KHz	550 kHz	350 - 600 kHz	450 KHz	3.8 MHz	3.8 MHz
Handswitching Handpiece	Yes, w/ power up/down	Yes	Yes, w/ power up/down	Yes, w/ power up/down	N/A	Optional
Footswitching Handpiece	Handswitching Handpiece works w/ optional footswitch	Handswitching Handpiece works w/ optional footswitch	Handswitching Handpiece works w/ optional footswitch	Optional - Purchase Handpiece & footswitch separately	Yes, Comes w/unit	Yes, Comes w/unit
Patient Plate	Optional	Optional	Yes - required for certain uses	Optional	Yes - required for all uses	Yes - required for all uses
Electrode Size	1/16" & 3/32"	3/32"	3/32"	3/32"	1/16" & 3/32"	1/16" & 3/32"
Warranty	2 yrs on unit	2 yrs on unit	2 yrs on unit	1 yr on unit	1 yr on unit	1 yr on unit