

## Step by step guide to the basics of taking dermatological pictures with the Nikon Coolpix 8400

### Mode

Indicator in the monitor: **S**

Mode dial up S → select 1/60 with mode dial, aperture is self-acting

### Image size

Indicator in the monitor: **2M**

MENU button → MY MENU → Image mode → Image size → 2M (1600x1200)

### Image sharpening

MENU button → MY MENU → Image sharpening → High

### White balance

Indicator in the monitor: **PRE**

MENU button → MY MENU → White balance → White bal. preset → Measure

### Focus mode

Indicator in the monitor: **Flower icon**

Button Focus mode → Flower icon

### Flash mode

Indicator in the monitor: 

Button flash mode → flash cancel

