

Sterilizer Comparisons

Manufacturer	Tuttnauer™	Midmark	Pelton Crane	SciCan	Harvey	ADEC
Model	EZ9™	M9	Delta XL8	n/a	n/a	n/a
Chamber size	9" x 18"	9" x 15"	8" x 14"	-	-	-
Estimated dry to dry cycle time (cold)	55	52	58	-	-	-
Active closed door drying	yes	no	yes	-	-	-
Warranty	2 yr. parts & labor	1 yr. parts & labor	3 yr. parts No labor	-	-	-
Model	EZ10™	M11	Delta XL10	Quantim	PV Dry	n/a
Chamber size	10" x 19"	11" x 18"	10" x 18"	10" x 13"	10" x 16"	-
Estimated dry to dry cycle time (cold)	55	60	58	57	44-54	-
Active closed door drying	yes	no	yes	yes	yes	-
Warranty	2 yr. parts & labor	1 yr. parts & labor	3 yr. parts No labor	2 yr. parts In house labor	1 yr. parts 90 day labor	-
Model	EZ10K™	n/a	n/a	n/a	n/a	Lisa
Chamber size	10" x 19"	-	-	-	-	10" x 14"
Estimated dry to dry cycle time (cold)	45	-	-	-	-	47- 54
Active closed door drying	yes	-	-	-	-	yes
Warranty	2 yr. parts & labor	-	-	-	-	1 yr. parts or 1,000 cycles

- Tuttnauer has the deepest chamber of any tabletop autoclave
- Cassette capacity 3 full plus 3 half standard (EZ10 and EZ10k)
- EZ autoclaves are the only autoclaves with a 2 year parts and labor warranty plus a 10 year chamber warranty