

MATERIAL SAFETY DATA SHEET

BLUE TISSUE STAIN

Delasco Dermatologic Lab & Supply, Inc.
608 13th Ave., Council Bluffs, IA 51501 1-712-323-3269
www.delasco.com questions@delasco.com
Date prepared: December 2010

The statements herein are offered for informational purposes only and are intended to be followed only by persons having related technical skills and at their own discretion and risk. Since conditions and manner of use are outside our control, we make no warranties, express or implied, and assume no liability in connection with any use of this information.

Product/Chemical Name: BLUE TISSUE STAIN

Formula: Mrs. Stewart's Double Concentrated Liquid Bluing

Physical Description: Blue Liquid with slightly pungent odor.
Soluble in water.

Special Hazards: None known

Unusual Fire/Explosion Hazards: This product presents no unusual fire or explosion hazards.

Reactivity: Stable. Avoid mixing with chlorites or phyochlorites. Avoid elevated temperatures

Hazardous Decomposition Products: None known

Spill or Leak Procedures

Action: Dispose of according to applicable local, state and federal regulations.

Disposal Method: To be performed in compliance with all current local, state and federal regulations.

Threshold Limit Value: Formaldehyde 2 ppm OSHA PEL; 1 ppm ACGIH TLV

Overexposure Effects: No hazards of ingestion or inhalation known.

Potential Carcinogen: None

First Aid Procedures

Skin: Wash with soap or mild detergent and large amounts of water until no evidence of chemical remains (15-20 minutes). Seek immediate medical attention.

Eyes: Check for and remove any contact lenses. Flush with cool water and consult physician

Ingestion: Check for and remove any dentures. No hazards of ingestion known

Inhalation: No hazards of inhalation known

Storage Recommendations: Store protected from light at room temperature.

Special Protection Information

No special protection needed.

Special Handling Precautions

N/A