

MATERIAL SAFETY DATA SHEET

DETACHOL

Delasco Dermatologic Lab & Supply, Inc.
608 13th Ave., Council Bluffs, IA 51501 1-712-323-3269
www.delasco.com questions@delasco.com
Date prepared: December 2010

The statements herein are offered for informational purposes only and are intended to be followed only by persons having related technical skills and at their own discretion and risk. Since conditions and manner of use are outside our control, we make no warranties, express or implied, and assume no liability in connection with any use of this information.

Product: DETACHOL (Adhesive Remover)

Chemical Name / Formula: Paraffins (C_nH_{2n+2}) and Cycloparaffins ($C_nH_{2(n+1-g)}$)

Physical Description: Reddish orange liquid

Boiling Point: Varies Specific gravity: 0.75-0.77 Water insoluble

Special Hazards: N/A

Unusual Fire/Explosion Hazards: Keep work areas free of hot metal surfaces and other sources of ignition.

Reactivity: This product is incompatible with strong oxidizing agents, strong acids or bases and selected amines.

Hazardous Decomposition Products: Thermal decomposition in the presence of air may yield carbon monoxide and/or carbon dioxide.

Spill or Leak Procedures

Action: Shut off and eliminate all ignition sources. Recover free product. Add sand or other suitable absorbant material to spill area. Keep product out of sewers.

Disposal Method: To be performed in compliance with all current local, state and federal regulations.

Threshold Limit Value: N/A

Overexposure Effects: High vapor concentration (more than 1000 ppm) are irritating to the eyes and the respiratory tract. May cause headaches and dizziness.

Possible carcinogen: No

First Aid Procedures

Skin: Remove contaminated clothing and shoes and wash skin thoroughly with soap and water (approx. 15-20 minutes) until no evidence of chemical remains.

Eye: Check for and remove any contact lenses. Flush with water for 15-20 minutes or until irritation subsides. If irritation persists call physician.

Ingestion: Check for and remove any dentures. If ingested do not induce vomiting. Call physician or poison control center immediately.

Inhalation: If overcome by vapor, remove from exposure and call physician. If breathing is irregular or has stopped, administer CPR. Get medical attention immediately.

Storage Recommendations: Store at room temperature.

Special Protection Information

Use supplied air respiratory protection in an enclosed area.

Special Handling Precautions

Chemical resistant gloves and splash resistant safety goggles are required.