

MATERIAL SAFETY DATA SHEET

DELASCO INSTRUMENT SPRAY LUBRICANT

Delasco Dermatologic Lab & Supply, Inc.
608 13th Ave., Council Bluffs, IA 51501 1-712-323-3269
www.delasco.com questions@delasco.com
Date prepared: September 2011

The statements herein are offered for informational purposes only and are intended to be followed only by persons having related technical skills and at their own discretion and risk. Since conditions and manner of use are outside our control, we make no warranties, express or implied, and assume no liability in connection with any use of this information.

Product/Chemical Name: DELASCO INSTRUMENT SPRAY LUBRICANT

Formula: Not available

Physical Description: White, milk-like liquid, slight oleic odor, Polymeric Lubricant Concentrate

Specific Gravity: 1.00 Water soluble

Special Hazards: Non-hazardous

Unusual Fire/Explosion Hazards: None

Reactivity: Stable

Hazardous Decomposition Products: None

Spill or Leak Procedures

Action: Add absorbant, sweep up and place in plastic drums for disposal.

Disposal Method: To be performed in compliance with all current local, state and federal regulations.

Threshold Limit Value: N/A

Overexposure Effects: Unknown

Fire Procedures

General Hazzard: Low hazzard, liquid can burn upon heating to temperatures above the flashpoint.

Flashpoint: 293°

Action: Use self contained breathing apparatus. Use water spray, carbon dioxide, dry chemical powder or apporpriate foam.

Potential Carcinogen: No

First Aid Procedures

Skin: Not harmful, however remove contaminated clothing before use.

Eye: Check for and remove any contact lenses. Flush with large amounts of water until no evidence of chemical remains (approx. 15-20 minutes). Seek medical attention if irritation occurs.

Ingestion: Not expected to cause serious harm.

Inhalation: Move to fresh air if irritation occurs. Seek medical attention.

Storage Recommendations: Store at room temperature.

Special Protection Information

Wear protective gloves when using the product.