

# MATERIAL SAFETY DATA SHEET

## MASTISOL

**Delasco Dermatologic Lab & Supply, Inc.**  
608 13th Ave., Council Bluffs, IA 51501 1-712-323-3269  
[www.delasco.com](http://www.delasco.com) [questions@delasco.com](mailto:questions@delasco.com)  
Date prepared: January 2008

The statements herein are offered for informational purposes only and are intended to be followed only by persons having related technical skills and at their own discretion and risk. Since conditions and manner of use are outside our control, we make no warranties, express or implied, and assume no liability in connection with any use of this information.

**Product/Chemical Name: MASTISOL (Surgical Adhesive)**

**Formula:**

**Physical Description:** Golden color liquid with aromatic odor.

Boiling Point: 190°F Specific Gravity: 0.85-0.9 Water soluble - partial

**Special Hazards:**

**Unusual Fire/Explosion Hazards:** Firefighters should wear self-contained breathing apparatus in the positive pressure mode with a full facepiece when there is a possibility of exposure to smoke fumes or hazardous decomposition products.

**Reactivity:** Contact with a wide range of oxidizing agents may cause violent reaction.

**Hazardous Decomposition Products:** Carbon dioxide can form on incomplete combustion.

**Spill or Leak Procedures**

**Action:** Use absorbent material like sand. Dispose absorbant material in a land-fill in accordance with local, state and Federal Regulations.

**Disposal Method:** To be performed in compliance with all current local, state & Federal Regulations.

**Threshold Limit Value:** N/A

**Overexposure Effects:** Headache, drowsiness and inability to concentrate. Ingestion of large quantities may cause depression of central nervous system, nausea and vomiting.

**Potential Carcinogen:** Unknown

**First Aid Procedures**

**Skin:** Flush exposed skin with water, wash with soap until no evidence of chemical remains.

**Eyes:** Check for and remove any contact lenses. Wash eyes immediately with large amounts of water or normal saline, until no evidence of chemical remains (approx. 15-20 minutes). Seek medical attention.

**Ingestion:** Check for and remove any dentures. Contact physician or poison control center immediately.

**Inhalation:** Remove to fresh air or well ventilated area.

**Special Protection Information**

Organic canister mask recommended. Safety goggles and gloves recommended.

**Special Handling Precautions**

N/A