

MATERIAL SAFETY DATA SHEET

RED TISSUE STAIN, MERCURY-FREE

Delasco Dermatologic Lab & Supply, Inc.
608 13th Ave., Council Bluffs, IA 51501 1-712-323-3269
www.delasco.com questions@delasco.com

Date prepared: September 2011

The statements herein are offered for informational purposes only and are intended to be followed only by persons having related technical skills and at their own discretion and risk. Since conditions and manner of use are outside our control, we make no warranties, express or implied, and assume no liability in connection with any use of this information.

Product/Chemical Name: RED TISSUE STAIN, MERCURY-FREE

Formula: Not available

Physical Description: Red liquid

Special Hazards: Use NIOSH approved respirator if excessive misting occurs.

Unusual Fire/Explosion Hazards: Dried films may burn.

Reactivity: Stable

Hazardous Decomposition Products: None

Spill or Leak Procedures

Action: Contain spill. Absorb with inert material (e.g., dry sand or earth), then place in chemical waste container. Avoid runoff into storm sewers and ditches which lead to waterways.

Disposal Method: To be performed in compliance with all current local, state and federal regulations.

Threshold Limit Value: N/A

Overexposure Effects: Eye contact may cause irritation. Prolonged or repeated skin contact may cause irritation.

Potential Carcinogen: N/A

First Aid Procedures

Skin: Wash with soap and water. Get medical attention if irritation develops or persists.

Eye: Check for and remove any contact lenses. Immediately flush eyes with plenty of water. Get medical attention if irritation persists.

Ingestion: Check for and remove any dentures. Do NOT induce vomiting. Give victim a glass of water or milk. Call a physician or poison control center immediately. Never give anything by mouth to an unconscious person.

Inhalation: Remove to fresh air. If not breathing, give artificial respiration. If breathing is difficult, give oxygen. Get immediate medical attention.

Storage Recommendations: Store in a cool, dry place.

Special Protection Information

Safety glasses with side shields and chemical resistant gloves are required.

Special Handling Precautions

Keep from freezing. Wash hands before eating. Remove contaminated clothing and wash before reuse.