

MATERIAL SAFETY DATA SHEET

RED TISSUE STAIN

Delasco Dermatologic Lab & Supply, Inc.
608 13th Ave., Council Bluffs, IA 51501 I-712-323-3269
www.delasco.com questions@delasco.com
Date prepared: January 2008

The statements herein are offered for informational purposes only and are intended to be followed only by persons having related technical skills and at their own discretion and risk. Since conditions and manner of use are outside our control, we make no warranties, express or implied, and assume no liability in connection with any use of this information.

Product/Chemical Name: RED TISSUE STAIN

Formula: MERBROMIN contains 27% mercury, and the formula is $C_{20}H_8HgNa_2O_6$

Physical Description: Red liquid

Special Hazards: Self-contained breathing apparatus is required

Unusual Fire/Explosion Hazards: Negligible

Reactivity: Stable

Hazardous Decomposition Products: None

Spill or Leak Procedures

Action: Take up with sand or other absorbent material and place into clean, dry containers for later disposal.

Disposal Method: To be performed in compliance with all current local state & Federal Regulations.

Threshold Limit Value: Not available

Overexposure Effects: Beware - can cause Mercury poisoning

Potential Carcinogen: No

First Aid Procedure:

Skin: Not available

Eyes: Check for and remove any contact lenses.

Ingestion: Check for and remove any dentures. May cause mercury poisoning

Inhalation: May cause mercury poisoning

Storage Recommendations: Store in a cool, dry place.

Special Protection Information

Eye goggles and rubber gloves are required. Wear eye shield.

Special Handling Precautions

Handle with care. Toxic causes mercury poisoning.