

MATERIAL SAFETY DATA SHEET

TOLUIDINE BLUE O

Delasco Dermatologic Lab & Supply, Inc.
608 13th Ave., Council Bluffs, IA 51501 1-712-323-3269
www.delasco.com questions@delasco.com
Date prepared: January 2008

The statements herein are offered for informational purposes only and are intended to be followed only by persons having related technical skills and at their own discretion and risk. Since conditions and manner of use are outside our control, we make no warranties, express or implied, and assume no liability in connection with any use of this information.

Product/Chemical Name: TOLUIDINE BLUE O

Formula: C₁₅H₁₆N₃SCL

Physical Description: Blue powder. Soluble in water.

Special Hazards: Self-contained breathing apparatus is required

Unusual Fire/Explosion Hazards: May emit very toxic fumes of CL, SO_x, NO_x if heated to decomposition

Reactivity: Stable

Hazardous Decomposition Products: CL, SO_x, NO_x and carbon oxides

Spill or Leak Procedures

Action: Wear protective gear, collect the spill and put into a proper container for disposal.

Disposal Method: To be performed in compliance with all current local, state & Federal Regulations.

Threshold Limit Value: Not available

Overexposure Effects: Unknown

Potential Carcinogen: Not available

First Aid Procedures

Skin: Remove contaminated clothing. Flush area with water for at least 15 minutes. Seek medical attention.

Eyes: Check for and remove any contact lenses. Flush with water for 15 minutes. Do not use an eye ointment. Seek medical attention.

Ingestion: Check for and remove any dentures. Do NOT induce vomiting. Seek medical attention immediately.

Inhalation: Remove from exposure area. Seek medical attention

Storage Recommendations: Store protected from light at room temperature.

Special Protection Information

Safety goggles and gloves required when handling.

Special Handling Precautions

Avoid contact with skin and eyes or breathing dust.