

MATERIAL SAFETY DATA SHEET

YELLOW TISSUE STAIN

Delasco Dermatologic Lab & Supply, Inc.
608 13th Ave., Council Bluffs, IA 51501 1-712-323-3269
www.delasco.com questions@delasco.com
Date prepared: December 2010

The statements herein are offered for informational purposes only and are intended to be followed only by persons having related technical skills and at their own discretion and risk. Since conditions and manner of use are outside our control, we make no warranties, express or implied, and assume no liability in connection with any use of this information.

Product/Chemical Name: YELLOW TISSUE STAIN

Formula: N/A

Physical Description: Yellow liquid which is soluble in water and has a mild alcohol odor.

Special Hazards: None under normal conditions

Unusual Fire/Explosion Hazards: N/A

Reactivity: Stable

Hazardous Decomposition Products: N/A

Spill or Leak Procedures

Action: For small spills: Wipe up and store wipers in underwriters approved containers.

Disposal Method: To be performed in compliance with all current local, state and Federal Regulations.

Threshold Limit Value: N/A

Overexposure Effects: Prolonged contact with skin or eyes can cause irritation.

Potential Carcinogen: No

First Aid Procedures

Skin: Wash affected area with soap and water until no evidence of chemical remains.

Eyes: Check for and remove any contact lenses. Wash eyes immediately with large amounts of water or normal saline, occasionally lifting upper and lower lids until no evidence of chemical remains (approx. 15-20 minutes). Seek medical assistance as required.

Ingestion: N/A

Inhalation: N/A

Storage Recommendations: Store protected from light at room temperature.

Special Protection Information

Local exhaust ventilation desirable. Wear protective gloves when handling.

Special Handling Precautions

Avoid contact with skin, eyes and clothing. Store in closed containers.